

# Lumenis Slit Lamps Specifications

	Lumenis 980 Parallel Optics	Lumenis 1000 Parallel Optics, Fully Integrated	Lumenis 990 Convergent Optics	SL 130 Zeiss Original
<b>Magnification</b>	5X	5X	2X / 3X / 5X	5X
<b>Compatible Devices</b>	Smart532™ Novus Spectra Vision One™ Array™ LaserLink™	Smart532™ Novus Spectra Vision One™	Selecta II Smart532™ Novus Spectra Vision One™	Smart532™ Novus Spectra Vision One™ Array™ LaserLink™
<b>LaserLink*</b>	LaserLink Z: 50-1000 µm spot size	LaserLink Z: 50-1000 µm spot size	LaserLink HS: 50-500 µm	LaserLink Z: 50-1000 µm spot size
<b>Eye Safety Filter (ESF)**</b>	Choice of automatic or fixed	Built-in automatic	Fixed, part of the LaserLink	Choice of automatic or fixed
<b>SureSpot™ Optics</b>	✓	✓	✓	✓
<b>Table</b>	Wheelchair accessible with adjustable height and locking castor wheels			
<b>Input Power</b>	Universal 100-240V			
<b>Micromanipulator</b>	Part of LaserLink Z	Part of slit lamp	None	Part of LaserLink Z

\* The LaserLink is built-in for Lumenis 1000 and is an optional addition to Lumenis 980, 990 and SL 130.

\*\* The eye safety filter is built-in for Lumenis 1000 and is a required addition to Lumenis 980 and SL 130 in order to perform laser treatments. For Lumenis 990, the ESF is built-in LaserLink HS.

## Lumenis Slit Lamp Portfolio



**Lumenis Inc.**  
1870 South Milestone Drive  
Salt Lake City, UT 84104 USA  
Tel +801-656-2300  
Fax +801-656-2429

**Lumenis (Germany) GmbH**  
Heinrich-Hertz-Str 3 D-63303  
Dreieich-Dreieichenhain  
GERMANY  
Tel +49 (0) 6103 8335 0

**AMERICAS**  
San Jose, CA, USA  
Tel +1 408 764 3000  
+1 877 586 3647  
Fax +1 408 764 3999

**EMEA**  
Dreieich Dreieichenhain,  
Germany  
Tel +49 6103 8335 0  
Fax +49 6103 8335 300  
Roma (RM), Italy  
Tel +39 06 90 75 230  
Fax +39 06 90 75 269  
Hertfordshire, UK  
Tel +44 20 8736 4110  
Fax +44 20 8736 4119

**JAPAN**  
Tokyo, Japan  
Tel +81 3 4431 8300  
Fax +81 3 4431 8301

**ASIA / PACIFIC**  
Beijing, China  
Tel +86 10 5737 6677  
Gurgaon, India  
Tel +91 124 422 07 95  
Kowloon, Hong Kong  
Tel +852 217 428 00  
Fax +852 272 251 51

ISO 13485 CERTIFIED



[www.lumenis.com/Ophthalmology](http://www.lumenis.com/Ophthalmology)

Lumenis® Certified Service | USA Toll-free 1-877-LUMENIS (1-877-586-3647)

© 2017 All Rights Reserved. The Lumenis Group of Companies. PB-2008216 Rev. A

**Lumenis®**  
Energy to Eye Care

**Lumenis®**  
Energy to Eye Care

Lumenis offers you a wide selection of advanced slit lamps.  
Choose the slit lamp that best fits your needs!



## Lumenis 980

Parallel optics (Zeiss style) diagnostic slit lamp

- 5X magnification changer
- Compatible with Lumenis photocoagulators and pattern scanner: Smart532™, Novus Spectra, Vision One™ and Array™ LaserLink™
- Table: wheelchair accessible with adjustable height and locking castor wheels
- Universal 100-240V input power

Optional additions:

- Choice of automatic or fixed eye safety filters
- LaserLink Z: 50-1000 μm spot size, containing SureSpot™ optics and micromanipulator with finger rest

## Lumenis 1000

Parallel optics (Zeiss style) fully integrated diagnostic slit lamp

- 5X magnification changer
- Compatible with Lumenis photocoagulators: Smart532™, Novus Spectra and Vision One™
- Table: wheelchair accessible with adjustable height and locking castor wheels
- Universal 100-240V input power
- Built-in micromanipulator
- Transferable (right or left) elbow rest
- Built-in three-wavelength automatic eye safety filter (532nm, 577nm, 659nm)
- LaserLink Z: 50-1000 μm spot size, containing SureSpot™ optics

## Lumenis 990

Convergent optics (Haag-Streit style) diagnostic slit lamp

- Available in 2X, 3X and 5X magnification changer
- Compatible with Lumenis portable SLT and photocoagulators: Selecta II, Smart532™, Novus Spectra and Vision One™
- Tonometer mount included
- Table: wheelchair accessible with adjustable height and locking castor wheels
- Universal 100-240V input power

Optional addition:

- LaserLink HS: 50-500 μm spot size, containing SureSpot™ optics and fixed eye safety filter

## SL 130

Original Zeiss, parallel optics, diagnostic slit lamp

- 5X magnification changer
- Compatible with Lumenis photocoagulators and pattern scanner: Smart532™, Novus Spectra, Vision One™ and Array™ LaserLink™
- Table: wheelchair accessible with adjustable height and locking castor wheels
- Universal 100-240V input power

Required additions:

- Choice of automatic or fixed eye safety filters
- LaserLink Z: 50-1000 μm spot size, containing SureSpot™ optics and micromanipulator with finger rest